


MEDI PLUS



Sterilization

Mesh Trays & Baskets

ABOUT US

We offer a full range of Mesh Trays and Baskets, Surgical and Dental Instruments specially tailored to the need of our valued clients innovative thinking, professional approach and years of experience has enabled us to put all the sectors of production together like a jigsaw. For more than three decades, highly skilled and experienced craftsmen engaged to produce and serve quality Surgical and Dental Instruments at competitive prices have acquired good reputation across the world enrolling big names of the field as their satisfied associates. Focusing on ultimate leadership and serving humanity



**Choose the Right
Medical Baskets
and Trays**



Product Catagory

- Sterile & Non-Sterile Single Use Surgical Instruments
- Surgical Reusable Instruments
- Sterile & Non-Sterile Dental Instruments (Reusable & Single Use)
- Mesh Trays & Baskets
- Hollowware Instruments
- Diagnostics Instruments
- Laboratory Instruments
- New Product Sterile & Non-Sterile Single Use Kits/Packs



PERFORATED / MESH TRAYS AND BASKETS

These Trays and baskets are designed for perfect washing and Sterilizing of surgical instruments. These modular steel trays are manufactured from electrolytically polished AISI 301 stainless steel. This means that they are compatible with all type of detergents, including strong alkalis.

Base of the trays is available in 3 different Mesh.

- Woven (Plain) Stainless Steel Mesh
- Welded Stainless Steel Mesh
- Woven Micro Stainless Steel Mesh

Item Code	Size in mm
SRR-PWB-401	120 x 250 x 30
SRR-PWB-402	120 x 250 x 50
SRR-PWB-403	120 x 250 x 60
SRR-PWB-404	120 x 250 x 70
SRR-PWB-405	120 x 250 x 100
SRR-PWB-406	240 x 250 x 30
SRR-PWB-407	240 x 250 x 50
SRR-PWB-408	240 x 250 x 60
SRR-PWB-409	240 x 250 x 70
SRR-PWB-410	240 x 250 x 100
SRR-PWB-411	255 x 250 x 30
SRR-PWB-412	255 x 250 x 50
SRR-PWB-413	255 x 250 x 60
SRR-PWB-414	255 x 250 x 70
SRR-PWB-415	255 x 250 x 100
SRR-PWB-416	300 x 250 x 30
SRR-PWB-417	300 x 250 x 50
SRR-PWB-418	300 x 250 x 60
SRR-PWB-419	300 x 250 x 70
SRR-PWB-420	300 x 250 x 100



Item Code	Size in mm
SRR-PWB-421	405 x 250 x 30
SRR-PWB-422	405 x 250 x 50
SRR-PWB-423	405 x 250 x 60
SRR-PWB-424	405 x 250 x 70
SRR-PWB-425	405 x 250 x 100
SRR-PWB-426	405 x 255 x 30
SRR-PWB-427	405 x 255 x 50
SRR-PWB-428	405 x 255 x 60
SRR-PWB-429	405 x 255 x 70
SRR-PWB-430	405 x 255 x 100
SRR-PWB-431	480 x 360 x 30
SRR-PWB-432	480 x 360 x 50
SRR-PWB-433	480 x 360 x 60
SRR-PWB-434	480 x 360 x 70
SRR-PWB-435	480 x 360 x 100
SRR-PWB-436	480 x 480 x 30
SRR-PWB-437	480 x 480 x 50
SRR-PWB-438	480 x 480 x 60
SRR-PWB-439	480 x 480 x 70
SRR-PWB-440	480 x 480 x 100
SRR-PWB-441	515 x 245 x 30
SRR-PWB-442	515 x 245 x 50
SRR-PWB-443	515 x 245 x 60
SRR-PWB-444	515 x 245 x 70
SRR-PWB-445	515 x 245 x 100
SRR-PWB-446	540 x 250 x 30
SRR-PWB-447	540 x 250 x 50
SRR-PWB-448	540 x 250 x 60
SRR-PWB-449	540 x 250 x 70
SRR-PWB-450	540 x 250 x 100

■ WIRE MESH TRAYS & BASKETS

These Trays and baskets are designed for perfect washing and Sterilizing of surgical instruments. These modular steel trays are manufactured from electrolytically polished AISI 301 stainless steel. This means that they are compatible with all type of detergents, including strong alkalis.

Base of the trays is available in 3 different Mesh.

- Woven (Plain) Stainless Steel Mesh
- Welded Stainless Steel Mesh
- Woven Micro Stainless Steel Mesh

Item Code	Size in mm
SRR-WMB-601	120 x 250 x 30
SRR-WMB-602	120 x 250 x 50
SRR-WMB-603	120 x 250 x 60
SRR-WMB-604	120 x 250 x 70
SRR-WMB-605	120 x 250 x 100
SRR-WMB-606	240 x 250 x 30
SRR-WMB-607	240 x 250 x 50
SRR-WMB-608	240 x 250 x 60
SRR-WMB-609	240 x 250 x 70
SRR-WMB-610	240 x 250 x 100
SRR-WMB-611	255 x 250 x 30
SRR-WMB-612	255 x 250 x 50
SRR-WMB-613	255 x 250 x 60
SRR-WMB-614	255 x 250 x 70
SRR-WMB-615	255 x 250 x 100
SRR-WMB-616	300 x 250 x 30
SRR-WMB-617	300 x 250 x 50
SRR-WMB-618	300 x 250 x 60
SRR-WMB-619	300 x 250 x 70
SRR-WMB-620	300 x 250 x 100
SRR-WMB-621	405 x 250 x 30
SRR-WMB-622	405 x 250 x 50
SRR-WMB-623	405 x 250 x 60
SRR-WMB-624	405 x 250 x 70
SRR-WMB-625	405 x 250 x 100
SRR-WMB-626	405 x 255 x 30
SRR-WMB-627	405 x 255 x 50
SRR-WMB-628	405 x 255 x 60
SRR-WMB-629	405 x 255 x 70
SRR-WMB-630	405 x 255 x 100



Wire mesh Trays and baskets are available with special capping (protection strip) on demand.



Item Code	Size in mm
SRR-WMB-631	480 X 360 X 30
SRR-WMB-632	480 X 360 X 50
SRR-WMB-633	480 X 360 X 60
SRR-WMB-634	480 X 360 X 70
SRR-WMB-635	480 X 360 X 100
SRR-WMB-636	480 X 480 X 30
SRR-WMB-637	480 X 480 X 50
SRR-WMB-638	480 X 480 X 60
SRR-WMB-639	480 X 480 X 70
SRR-WMB-640	480 X 480 X 100
SRR-WMB-641	515 X 245 X 30
SRR-WMB-642	515 X 245 X 50
SRR-WMB-643	515 X 245 X 60
SRR-WMB-644	515 X 245 X 70
SRR-WMB-645	515 X 245 X 100
SRR-WMB-646	540 X 250 X 30
SRR-WMB-647	540 X 250 X 50
SRR-WMB-648	540 X 250 X 60
SRR-WMB-649	540 X 250 X 70
SRR-WMB-650	540 X 250 X 100

■ LID PERFORATED / MESH TRAYS & BASKETS

Separate Lids are also available.



Item Code	Size in mm
SRR-LID-501	120 x 250
SRR-LID-502	240 X 250
SRR-LID-503	255 X 245
SRR-LID-504	300 X 300
SRR-LID-505	405 X 250



Item Code	Size in mm
SRR-LID-501	405 x 255
SRR-LID-502	480 X 360
SRR-LID-503	480 X 480
SRR-LID-504	515 X 245
SRR-LID-505	540 X 250

■ LID WIRE MESH BASKETS

Separate Lids are also available.




Item Code	Size in mm
SRR-LID-501	120 x 250
SRR-LID-502	240 X 250
SRR-LID-503	255 X 245
SRR-LID-504	300 X 300
SRR-LID-505	405 X 250



Item Code	Size in mm
SRR-LID-506	405 x 255
SRR-LID-507	480 X 360
SRR-LID-508	480 X 480
SRR-LID-509	515 X 245
SRR-LID-510	540 X 250

■ MICRO / FINE MESH BASKETS

Mico Mesh Baskets are an ideal way of insuring small item are not lost or misplaced during the washing and sterilization cycle. Micromesh baskets are manufactured from electrolytically polished AISI304 Woven Stainless Stell Micro Mesh



Item Code	Size in mm
SRR-MMB-701	40 x 40 x 20
SRR-MMB-701	70 x 70 x 30
SRR-MMB-701	80 x 40 x 20
SRR-MMB-701	80 x 40 x 30
SRR-MMB-701	80 x 80 x 34
SRR-MMB-701	80 x 80 x 40
SRR-MMB-701	90 x 90 x 50
SRR-MMB-701	105 x 70 x 25
SRR-MMB-701	150 x 50 x 35
SRR-MMB-710	160 x 40 x 20
SRR-MMB-711	180 x 130 x 20
SRR-MMB-712	190 x 85 x 100
SRR-MMB-713	200 x 100 x 35
SRR-MMB-714	200 x 110 x 40
SRR-MMB-715	210 x 170 x 20
SRR-MMB-716	220 x 140 x 50
SRR-MMB-717	225 x 110 x 50
SRR-MMB-718	235 x 120 x 60
SRR-MMB-719	235 x 155 x 40
SRR-MMB-720	240 x 90 x 40
SRR-MMB-721	240 x 120 x 35
SRR-MMB-722	250 x 200 x 20
SRR-MMB-723	270 x 135 x 40
SRR-MMB-724	275 x 175 x 35
SRR-MMB-725	275 x 200 x 35
SRR-MMB-726	280 x 145 x 50
SRR-MMB-727	300 x 180 x 20
SRR-MMB-728	400 x 250 x 50

■ MICRO / FINE MESH BASKETS WITH CAPPING



Wire mesh Trays and baskets are available with special capping (protection strip) on demand.



Item Code	Size in mm
SRR-MMB-729	40 x 40 x 20
SRR-MMB-730	70 x 70 x 30
SRR-MMB-731	80 x 40 x 20
SRR-MMB-732	80 x 80 x 40
SRR-MMB-733	90 x 90 x 50
SRR-MMB-734	160 x 40 x 20
SRR-MMB-735	200 x 40 x 20
SRR-MMB-736	200 x 110 x 40
SRR-MMB-737	225 x 100 x 50
SRR-MMB-738	235 x 120 x 60
SRR-MMB-739	270 x 135 x 40
SRR-MMB-739	280 x 145 x 50

■ FINE MESH BASKET WITH SEPARATE LID




Item Code	Size in mm
SRR-FMB-801	120 X 155 X 40
SRR-FMB-802	120 X 155 X 50
SRR-FMB-803	235 X 155 X 40
SRR-FMB-804	235 X 155 X 50
SRR-FMB-805	375 X 260 X 40
SRR-FMB-806	380 X 265 X 45
SRR-FMB-807	380 X 265 X 55

■ FINE MESH BASKET WITH REMOVABLE LOCKEABLE LID

Item Code	Size in mm
SRR-FMB-808	120 X 155 X 40
SRR-FMB-809	120 X 155 X 50
SRR-FMB-810	235 X 155 X 40
SRR-FMB-811	235 X 155 X 50
SRR-FMB-812	375 X 260 X 40
SRR-FMB-813	380 X 265 X 45
SRR-FMB-814	380 X 265 X 55



■ FINE MESH BASKET WITH CATCH LOCK



Item Code	Size in mm
SRR-FMB-815	235 X 80 X 30
SRR-FMB-816	280 X 80 X 30

■ ENDSCOPE / ARTHROSCOPY BASKET



Item Code
SRR-ESB-903
SRR-ESB-904

Size in mm
290 x 80 x 52
290 x 80 x 52
(Double Scope)


Item Code
SRR-ESB-901
SRR-ESB-902

Size in mm
180 x 80 x 52
180 x 80 x 52
(Double Scope)



Item Code
SRR-ESB-905
SRR-ESB-906

Size in mm
290 x 80 x 52
290 x 80 x 52
(Double Scope)



Item Code
SRR-ESB-907
SRR-ESB-908
SRR-ESB-909

Size in mm
600 x 80 x 52
670 x 80 x 52
700 x 80 x 52


■ ENDOSCOPE / ARTHROSCOPY BASKET



Item Code	Size in mm
SRR-ECB-910	290 x 160 x 52
SRR-ECB-911	460 x 160 x 52



Item Code	Size in mm
SRR-ECB-912	290 x 160 x 52
SRR-ECB-913	460 x 160 x 52



Item Code	Size in mm
SRR-ECB-914	460 x 160 x 52

■ STERILE GOODS BASKETS



Item Code	Size in mm
SRR-EGB-1001	575 x 280 x 135
SRR-EGB-1002	575 x 280 x 265

Item Code	Size in mm
SRR-EGB-1003	575 x 280 x 265 (Front Access)



■ BASKET SPRI

Double frame with folding handle



Item Code	Size in mm
SRR-EGB-1004	575 x 280 x 265
SRR-EGB-1005	585 x 395 x 195
SRR-EGB-1006	600 x 280 x 138
SRR-EGB-1007	600 x 395 x 260
SRR-EGB-1008	600 x 300 x 97
SRR-EGB-1009	600 x 300 x 138
SRR-EGB-1010	600 x 300 x 260

Single Frame

Item Code	Size in mm
SRR-EGB-1011	600 x 400 x 100
SRR-EGB-1012	600 x 400 x 130
SRR-EGB-1013	600 x 400 x 150
SRR-EGB-1014	600 x 400 x 180
SRR-EGB-1015	600 x 400 x 200



■ Toost Rack

Item Code
SRR-TRS-001
SRR-TRL-002



New Addition

■ TRINKET BOX PERFORATED SIDES WITH WIRE MESH LID

Item Code	Size in mm
SRR-TPB-2001	90 x 60 x 25
SRR-TPB-2002	105 x 70 x 25
SRR-TPB-2003	120 x 80 x 25
SRR-TPB-2004	150 x 100 x 25
SRR-TPB-2005	151 x 130 x 25



■ TRINKET BOX PERFORATED SIDES





Item Code	Size in mm
SRR-TPB-2006	90 x 60 x 25
SRR-TPB-2007	105 x 70 x 25
SRR-TPB-2008	120 x 80 x 25
SRR-TPB-2009	150 x 100 x 25
SRR-TPB-2010	151 x 130 x 25

■ TRINKET ROUND BASKETS

Item Code	Size in mm
SRR-TPB-2011	80 X 35
SRR-TPB-2012	90 X 35



■ INSTRUMENTS CASSETTS



Item Code	Size in mm
SRR-ICB-3001	220 X 125 X 40
SRR-ICB-3002	250 X 220 X 40

■ INSTRUMENTS CASSETTS



Item Code	Size in mm
SRR-ICB-3003	240 X 90 X 40
SRR-ICB-3004	235 X 80 X 30
SRR-ICB-3005	235 X 155 X 40



Item Code	Size in mm
SRR-ICB-3006	205 X 100 X 30
SRR-ICB-3007	205 X 125 X 30
SRR-ICB-3008	205 X 152 X 30
SRR-ICB-3009	205 X 230 X 30
SRR-ICB-3010	205 X 280 X 30



Item Code	Size in mm
SRR-ICT-3001	180 X 60 X 22
SRR-ICT-3002	180 X 125 X 22



Item Code	Size in mm
SRR-ICT-3003	184 X 74 X 34
SRR-ICT-3004	184 X 145 X 34
SRR-ICT-3005	205 X 143 X 34



Item Code	Size in mm
SRR-ICT-3006	280 X 183 X 32
SRR-ICT-3007	280 X 205 X 32

■ INSTRUMENTS CASSETTES



Item Code
SRR-ICT-3008

Size in mm
280 x 183 x 32

■ MICRO CASSETTES



Item Code	Size in mm
SRR-MCT-3009	185 x 82 x 18
SRR-MCT-3010	185 x 140 x 18



Item Code	Size in mm
SRR-MCT-3011	40 x 55 x 20
SRR-MCT-3012	74 x 45 x 20
SRR-MCT-3013	74 x 59 x 24

■ DENTAL BOX

Item Code	Size in mm
SRR-EDB-3014	250 x 160 x 40



■ ENDODONTIC BOX

Item Code	Size in mm
SRR-EDB-3015	105 x 45 x 50



■ IMPLANTOLOGY



Item Code
SRR-IMP-4001

Size in mm
180 X 148 X 50



Item Code
SRR-IMP-4002

Size in mm
80 X 55 X 40



Item Code
SRR-IMP-4003

Size in mm
150 X 85 X 50



Item Code
SRR-IMP-4004

Size in mm
100 X 90 X 50



Item Code
SRR-IMP-4005

Size in mm
150 X 85 X 50



Item Code
SRR-IMP-4006

Size in mm
80 X 55 X 40



Item Code
SRR-IMP-4007

Size in mm
100 X 90 X 50




Item Code
SRR-IMP-4008

Size in mm
80 X 55 X 40

■ ACCESSORIES

Sieves

Item Code	Size in mm
SRR-ACC-5001	140
SRR-ACC-5002	200




■ NEEDLE CASES



Item Code	Size in mm
SRR-ACC-5003	140 X 90 X 15



Item Code	Size in mm
SRR-ACC-5004	65 X 15



Item Code	Size in mm
SRR-ACC-5005	50 X 30 X 8
SRR-ACC-5006	65 X 30 X 8
SRR-ACC-5007	80 X 30 X 8

■ TEST TUBE STAND



Item Code	Size in mm
SRR-ACC-5008	12 Tubes nos.


■ ACCESSORIES

Fixing Calamps





Item Code	Size in mm
SRR-ACC-5009	60/220 X 220 X 39

Instruments Stringer




Item Code	Size in mm
SRR-ACC-5010	4" X 2-1/2"
SRR-ACC-5011	6" X 2-1/2"
SRR-ACC-5012	8" X 2-1/2"
SRR-ACC-5013	10" X 2-1/2"
SRR-ACC-5014	12" X 2-1/2"

Silicon Mats



Item Code	Size in mm
SRR-ACC-5015	515 X 240 X 15
SRR-ACC-5016	545 X 565 X 19

STAINLESS STEEL TAGS



Item Code	Size in mm
SRR-ACC-5020	50 X 25
SRR-ACC-5021	75 X 30



Item Code	Size in mm
SRR-ACC-5022	90 X 50
SRR-ACC-5023	90 X 35

NOMENCLATURE PLATE



Item Code	Size in mm
SRR-ACC-5017	35 X 35
SRR-ACC-5018	90 X 30
SRR-ACC-5019	90 X 40

■ HOSPITAL WARE

Perforated sheet Trays



Item Code	Size in mm
SRR-HPW-6001	255 X 255 X 50
SRR-HPW-6002	255 X 255 X 75
SRR-HPW-6003	255 X 255 X 100
SRR-HPW-6004	255 X 255 X 50
SRR-HPW-6005	480 X 255 X 75
SRR-HPW-6006	480 X 255 X 100
SRR-HPW-6007	510 X 255 X 50
SRR-HPW-6008	510 X 255 X 75
SRR-HPW-6009	510 X 255 X 100

Instruments Trays



Item Code	Size in mm
SRR-HPW-6010	200 X 100 X 25
SRR-HPW-6011	200 X 150 X 40
SRR-HPW-6012	230 X 150 X 50
SRR-HPW-6013	250 X 140 X 50
SRR-HPW-6014	250 X 200 X 50
SRR-HPW-6015	300 X 200 X 50
SRR-HPW-6016	300 X 250 X 50
SRR-HPW-6017	350 X 250 X 50
SRR-HPW-6018	350 X 300 X 50
SRR-HPW-6019	400 X 300 X 50
SRR-HPW-6020	450 X 300 X 50

Catheter Trays



Item Code	Size in mm
SRR-HPW-6021	230 X 80 X 50
SRR-HPW-6022	250 X 100 X 40
SRR-HPW-6023	300 X 100 X 40
SRR-HPW-6024	400 X 100 X 50

Instruments Box



Item Code	Size in mm
SRR-HPW-6025	150 X 80 X 40
SRR-HPW-6026	170 X 70 X 30
SRR-HPW-6027	180 X 80 X 40
SRR-HPW-6028	200 X 100 X 40
SRR-HPW-6029	230 X 120 X 50
SRR-HPW-6030	250 X 125 X 50
SRR-HPW-6031	300 X 150 X 50
SRR-HPW-6032	300 X 200 X 50
SRR-HPW-6033	350 X 175 X 50
SRR-HPW-6034	350 X 180 X 60
SRR-HPW-6035	400 X 200 X 90
SRR-HPW-6036	500 X 200 X 120

■ HOSPITAL WARE

Scaller Trays



Item Code	Size in mm
SRR-HPW-6038	200 x 100
SRR-HPW-6039	205 x 115
SRR-HPW-6040	210 x 150
SRR-HPW-6041	220 x 130
SRR-HPW-6042	225 x 100

Kidney Trays



Item Code	Size in mm
SRR-HPW-6043	150
SRR-HPW-6044	180
SRR-HPW-6045	200
SRR-HPW-6046	200
SRR-HPW-6047	220
SRR-HPW-6048	240
SRR-HPW-6049	250
SRR-HPW-6050	280



Item Code	Lid Size in mm
SRR-HPW-6051	150
SRR-HPW-6052	180
SRR-HPW-6053	200
SRR-HPW-6054	250
SRR-HPW-6055	280

Irrigators



Item Code	Size in Ltr
SRR-HPW-6056	1.0
SRR-HPW-6057	1.5
SRR-HPW-6058	2.0

Measuring Jugs



Item Code	Size in ml
SRR-HPW-6059	250
SRR-HPW-6060	500
SRR-HPW-6061	1000

■ HOSPITAL WARE

Dressing Drum



Item Code	Size in mm
SRR-HPW-6062	80 x 80
SRR-HPW-6063	120 x 120
SRR-HPW-6064	130 x 290
SRR-HPW-6065	150 x 140
SRR-HPW-6066	160 x 120
SRR-HPW-6067	180 x 180
SRR-HPW-6068	200 x 190
SRR-HPW-6069	240 x 240
SRR-HPW-6070	290 x 240
SRR-HPW-6071	340 x 240
SRR-HPW-6072	390 x 200
SRR-HPW-6073	390 x 290
SRR-HPW-6074	390 x 390

Forceps Jar



Item Code	Size in mm
SRR-HPW-6075	30 x 80
SRR-HPW-6076	180 x 50
SRR-HPW-6077	200 x 65
SRR-HPW-6078	200 x 100

Waste Collector



Item Code	Size in mm
SRR-HPW-6079	Standard

Lotion Bowls



Item Code	Size in mm
SRR-HPW-6080	80 x 30
SRR-HPW-6081	100 x 60
SRR-HPW-6082	160 x 60
SRR-HPW-6083	180 x 70
SRR-HPW-6084	230 x 80
SRR-HPW-6085	240 x 100
SRR-HPW-6086	250 x 100
SRR-HPW-6087	300 x 100
SRR-HPW-6088	350 x 125
SRR-HPW-6089	400 x 150

Gallipots



Item Code	Size in oz.
SRR-HPW-6099	2
SRR-HPW-6100	3
SRR-HPW-6101	4
SRR-HPW-6102	6
	8

Lodine Cups



Item Code	Size in mm
SRR-HPW-6090	80 x 40
SRR-HPW-6091	100 x 50
SRR-HPW-6092	160 x 50
SRR-HPW-6093	180 x 50
SRR-HPW-6094	230 x 75
SRR-HPW-6095	240 x 100
SRR-HPW-6096	250 x 100
SRR-HPW-6097	300 x 100
SRR-HPW-6098	350 x 115



Item Code	Size in oz.
SRR-HPW-6104	2
SRR-HPW-6105	3
SRR-HPW-6106	4
SRR-HPW-6107	6
	8

■ HOSPITAL WARE

Dressing Jar



Item Code	Size in mm
SRR-HPW-6109	80 x 80
SRR-HPW-6110	80 x 90
SRR-HPW-6111	100 x 80
SRR-HPW-6112	100 x 100
SRR-HPW-6113	120 x 120
SRR-HPW-6114	125 x 100
SRR-HPW-6115	150 x 150

Bed pan

Perfection



Item Code	Size in mm
SRR-HPW-6116	Adult
SRR-HPW-6117	Child

Slipper



Item Code	Size in mm
SRR-HPW-6118	Adult
SRR-HPW-6119	Child

Urinals

Male Urinal



Item Code	Size in mm
SRR-HPW-6120	Standard



Item Code	Size in mm
SRR-HPW-6121	Standard

Urinal Jug



Item Code	Size in mm
SRR-HPW-6122	Standard

Single Use Surgical Instruments



Dental Instruments



Diagnostic Instruments



Otoscopes



General Surgery Instruments



TC Instruments





MediPluss

Airport Road Kotli Behram,

Sialkot 51310-Pakistan

P +92 302 8712217 | E inquiry@medipluss.com

www.medipluss.com

